ACE DA-7100

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Important Notes

On the Calibration

Each breathalyzer which you buy from ACE Instruments can also be calibrated by ACE Instruments. Of course, we also calibrate devices which were bought from other contractors, as long as these devices are identical in construction to the devices sold by ACE Instruments.

You can see the due date of the next calibration on the calibration sticker placed on your breathalyzer. If the sticker is missing, the calibration is to be done 6 months after the purchase date: e.g. purchase date: 13.05. = next calibration on the 13.11. and the next on the 13.05. of the following year.

For a calibration service you can send the device to the service address specified at www.alkomat.net.

Why is a calibration necessary?

The calibration is not mandatory, but the accuracy of the sensor installed in the breathalyzer will decrease steadily without a calibration and additionally guarantee expires, if you don't pay attention to the calibration ranges.

Which factors influence the measurement accuracy

- The more often the device is used the more the zero point is influenced.
- The higher the alcohol concentration is, the more the zero point is relocated.
- Devices which are used rarely have a decreased sensibility, too.

ACE Instruments recommends using the breathalyzer at least once (1x) a week independently of having drunken alcohol or not. To maintain the measurement accuracy, the breathalyzer should also be calibrated periodically.

Please also be aware of the several factors of influence on the measuring result

- Alcohol amount
- Period of the consumption of alcohol
- Gender
- Age
- Weight
- Height
- Drinking habits
- Health state
- Vitality of the liver
- Size of the lung
- Admission of food during the consumption of alcohol
- Influence of medicine

Important notes on test execution and self-assessment

After the last sip of alcohol you should wait at least 15 minutes, as there may still be residual alcohol in the oral cavity, which influences the measurement result. Also smoking should be avoided before a test, as cigarette smoke harms the sensor.

For each test you should use a new mouthpiece – not only for hygienic reasons, but also because an already used mouthpiece can influence the next measurement result. The mouthpiece could contain residual alcohol, which would affect the new test result when the same mouthpiece is used again (this for example could show an alcoholization of a person who hasn't drunken any alcohol).

If "FLo" appears on the display, the given breath sample was insufficient or ended too soon. Blow into the mouthpiece for about 4 seconds constantly (without a break) – preferably as if you were whistling without a sound.

If the same person performs 2 tests consecutively, then wait at least 3 minutes between the two tests, as a too short waiting time may cause different measurement results. The air exhaled from the lungs is, after a test, no longer enriched with the same alcohol content as in the previous test. Therefore, an immediate second test will indicate a lower test result.

All measured values determined with this breathalyzer are always indicative and not legally binding.

Why the estimated alcoholization and the measurement result in most cases cannot match

An alcohol test is affected by a lot of different factors which also affect the test result. Beside the gender, age, body size and body weight, also the body fat percentage and the amount of body fluid have a big influence.

In addition, the general physical and mental state or the current daily constitution have to be taken into consideration. Furthermore, the meals eaten before a breath test (type of content and level of the stomach) and the intake of medication may influence the measurement.

Other influencing factors are, of course, the type of drink (percentage of alcohol content in the beverage and its carbon dioxide and sugar content), the drinking rate and the individual subject's duration of alcohol build-up or degradation in the body.

Some people quickly have a high alcohol content in the blood, with other people it takes longer until an alcoholization is detected. About 40-50% of the consumed alcohol is in the body within 15 minutes, after 30 minutes it is 65% and after 70 minutes the previously consumed alcohol is 100% absorbed in the body. However, it should also be noted that some people – partly due to genetic factors – quickly reduce the alcohol content in the blood, whereas others still have an increased alcohol content in the blood after a longer period of time.

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Operating Manual - EN

Author's Liability

All information herein was created conscientiously and under the assumption of correctness. Hence, the author cannot be held liable for the use of any kind whatsoever and/or any ensuing consequences. Furthermore, the author is not liable for any costs or expenses for damage that result from any information in this document. Any information is subject to change without prior notice.

Introduction

Thank you for purchasing our ACE DA-7100 Breathalyzer! It measures the alcohol concentration in human breath using a fuel cell sensor and calculates the blood alcohol level. The test results are very reliable and have a tolerance of +/- 0.05‰ (per mill) with reference to 1.00‰ at 25°C ambient temperature. The measurement results are displayed as a 3-digit number with 2 decimals and are interpretable 1:1 (for example the display shows the value 0.36 = 0.36 per mill). With the measurement unit % BAC the test result is displayed with 4 digits and 3 decimals.

The fuel cell sensor inside the Breathalyzer can be calibrated. To maintain the measurement accuracy of this device and the sensor warranty, it should be calibrated regularly. For a calibration service you can send the device to the service address specified at www.alkomat.net.

The device must only be used pursuant to the instructions mentioned in this manual. All measured values determined by this breathalyzer are always a guide only and are not legally binding. Breath tests should be carried out using clip-on mouthpieces (one new mouthpiece should be used per test for reasons of hygiene and to avoid falsified readings).

Note to country-specific measurement units

Customers in the EU get the breathalyzer with ‰ (per mill) as default setting, in the UK the default setting is mg/L and in the USA it is % BAC. In the chapter „Technical data“ you will find information about the accuracy concerning the different measurement units. Please check e.g. the website of the department of transportation, the website of the Ministry of Transport or of automotive clubs to get informed about the current alcohol limits in your country.

Scope of delivery

- 1 Digital Fuel Cell ACE DA-7100 Breathalyzer
- 2 AAA Batteries (1.5 V, Alkaline)
- 6 Mouthpieces
- 1 Carrying Pouch
- 1 Quick Start Guide (German & English)
- 1 Operating Manual (Multilingual)

Label of individual components

1. Mouthpiece and Mouthpiece Slot
2. Power Button
3. Digital LED Display (3 digits with unit ‰ and mg/L and 4 digits with unit % BAC)
4. Battery Compartment Cover (at the back)
5. Air Outlet (always keep this part unblocked)

This way you receive results that are as accurate as possible

- A) Wait at least 15 minutes after your last sip of alcohol. Do not eat, drink or smoke during that time.
- B) Do not open the casing.
- C) Do not keep the breathalyzer in cold, hot, moist or badly ventilated places.
- D) Keep the breathalyzer out of the reach of children.
- E) Conduct three inlet tests after long storage before conducting the actual alcohol test.
- F) Do not drop the breathalyzer.
- G) Change the battery if it is depleted. Diminished battery capacity leads to incorrect test results.
- H) Do not conduct tests in loud (sound waves disrupt the sensor), drafty or poorly ventilated places.
- I) After a test result of 1.000‰ (per mill = default setting) or higher the sensor needs rest time. In this case wait at least three minutes before the next test.
- J) Do not conduct tests in very cold, very hot or an unventilated, smoky ambience.
- K) Store the breathalyzer in a way to avoid penetration by dust or liquids.
- L) Blowing too hard may cause damage to the sensor as saliva flows into it.

Preparation for the test

1. If you use the device for the first time you have to remove the small paper inside the battery compartment which prevents the battery from discharging. Open the battery compartment cover at the back of the device (push it downwards), pull out the paper and close the battery compartment cover again. If you exchange old batteries against new ones please mind the battery polarity (+/- pole).
2. Insert a mouthpiece into the slot on the side of the breathalyzer. Ensure whether it securely fits or not. For sanitary reasons it is highly recommended to use an extra mouthpiece for each test.

Conducting the measurement

1. Push the power button, then you hear one beep sound and shortly will see the number of all tests performed with the device.
2. After pushing the power button a countdown starts (shown on the display) and the device warms up.
3. At the end of the countdown you see "0" on the display and additionally the device beeps twice. As soon as you can see "0", start blowing into the mouthpiece. Tip: Blow like you would whistle (without a break). The device will start beeping when you start blowing.
4. Keep blowing steadily until the beeping stops and you hear a double beep sound. If the breath sample was not sufficient you will see „FLo“ on the display. Push the start button and start a new test.
5. After a successful test the result is displayed after a few seconds. The test result is shown as 3 digit readout of your blood alcohol level in ‰ (per mill = EU setting) OR mg/L (= UK setting) OR % BAC (= US setting). If the test result is higher than 0.50‰ or 0.25mg/L or 0.05% BAC then you can hear a warn signal (beep sound) for about 10 seconds.

Please note: The device puts itself off automatically approx. 15 seconds after showing the test result or 10 seconds after not using it – additionally you hear a beep sound and you will see "Out" on the display.

The accuracy of this device is of +/- 0.05‰ (per mill = default setting) with reference to 1.00‰. For devices with UK measurement unit the tolerance is +/- 0.025mg/L with reference to 0.50mg/L, for devices with US measurement unit the tolerance is +/- 0.005% BAC with reference to 0.100% BAC.

ACE Instruments recommends

- ALWAYS USE A NEW MOUTHPIECE FOR EACH PERSON AND TEST THAT IS CONDUCTED. ACE INSTRUMENTS ASSUMES NO LIABILITY FOR ANY DAMAGES WHICH ARISE WHEN SEVERAL PEOPLE USE THE SAME MOUTHPIECE FOR MULTIPLE TESTS.
- ONLY USE THE PROVIDED SINGLE PACKED AND HYGIENIC MOUTHPIECES - MORE ARE AVAILABLE AT WWW.ALKOMAT.NET.
- DURING EACH TEST THE BREATHALYZER'S INTERNAL AND EXTERNAL PARTS COME INTO CONTACT WITH HUMAN SECRETION (SALIVA, BLOOD, SWEAT, ETC.)! USE DISPOSABLE GLOVES WHEN PULLING OUT A MOUTHPIECE USED BY ANOTHER PERSON! YOU CAN BUY APPROPRIATE GLOVES E.G. AT ACE TECHNIK (WWW.ACE-TECHNIK.DE).

IMPORTANT NOTICE

WHEN YOU HAVE EATEN, DRUNKEN ALCOHOL OR SMOKED A CIGARETTE, WAIT AT LEAST 15 MINUTES BEFORE STARTING THE MEASUREMENT. ANY TEST CONDUCTED PRIOR TO THE EXPIRATION OF THE RECOMMENDED WAITING PERIOD WILL RESULT IN INACCURATE RESULTS.

IF THE TEST RESULT IS BELOW 0.10‰ (PER MILL) OR 0.050 mg/L (UK setting) OR 0.010% BAC (US setting) THEN THE DEVICE SHOWS 0.00 (‰ OR mg/L OR % BAC – depending on the setting) ON THE DISPLAY. THIS HAPPENS BECAUSE AT A LOW RANGE OF ALCOHOL CONCENTRATION IT IS POSSIBLE TO MEASURE OTHER MATERIALS IN YOUR MOUTH.

This breathalyzer is a high quality and sensitive measurement device. To maintain the accuracy of the sensor, it should regularly be used. If the device is not used for an extended period of time, the sensor can dry out. ACE Instruments recommends to use the breathalyzer at least 2 times per month (even if you have not been drinking alcohol).

Due to changes in the product, the procedures stated can change. This can take place without prior notice or adaptation of the operating manual and is no cause for complaint, return or price reduction. The user is responsible for recognizing discrepancies between the operating manual and the actual handling of the breathalyzer and must react accordingly.

Only the usage of original ACE Instruments mouthpieces

- Guarantees test result accuracy
- Guarantees 100% hygiene
- Guarantees full warranty of this device

ACE INSTRUMENTS ADVISES

EVEN IF YOU GET A RESULT MEANING „SAFE RANGE“, YOU SHOULD NEVER DRIVE AFTER CONSUMING ALCOHOL OR DO ANY OTHER ACTIVITIES THAT ARE SEVERELY IMPAIRED BY ALCOHOL CONSUMPTION AND MIGHT ENDANGER YOUR HEALTH AND SAFETY OR THE HEALTH AND SAFETY OF OTHERS.

Display indications

Display	Meaning	Action
	Full battery	None required.
	Depleted battery: Warning signal when battery is low. The device is operational, but to ensure correct results, you should insert new batteries.	Insert new batteries.
BAT	When 'BAT' is displayed, the device is not operational. Change batteries immediately.	Insert new batteries immediately.
FLo	This is displayed when you have blown wrong or too weakly into the mouthpiece.	Press the start button again to repeat the test.
Out	When you have not blown into the mouthpiece within 30 seconds, this is displayed.	Press the start button to reactivate the countdown.
Out C	This is displayed when the device was started outside of the approved operating temperature (0 - 40° C).	Repeat the test as soon as the correct operating temperature is reached.

Calibration

To maintain the accuracy of the breathalyzer, ACE Instruments recommends a calibration after 500 tests or approx. every 12 months (depending which incident happens first). For a calibration service you can send the device to ACE Instruments (see www.alkomat.net). The calibration is not mandatory, but the accuracy of the breathalyzer will decrease steadily without it.

We recommend sending it as a parcel, so that in cases of loss an inquiry or proof of delivery is possible.

Due to changes in the product, the procedures stated can change. This can take place without prior notice or adaptation of the operating manual and is no cause for complaint, return or price reduction. The user is responsible for recognizing discrepancies between the operating manual and the actual handling of this breathalyzer and must react accordingly.

Battery notices

The device associated with this operating manual contains two AAA Alkaline batteries (1.5 V) in delivery condition. These batteries are not made for recharging. The use of rechargeable batteries is not possible with this device and can damage it. If the included batteries are depleted, they must not be thrown into the household waste. Used batteries may contain harmful substances that can cause environmental or health damage. Please return all depleted batteries to electrical appliance stores, who collect them for free.

Battery exchange

To replace the batteries, push the battery cover on the back of the unit downwards and replace the discharged batteries with two new ones. Be aware of the polarity when inserting the batteries (+/- pole). Put the battery compartment cover on the rear panel and press it gently on the device until it clicks. In the EU the return of used batteries is mandatory by law. All depleted batteries will be recycled. This way valuable resources like iron, zinc and nickel can be retrieved and are reusable then. Battery recycling is the easiest lesson in environmental protection – please support it!

The shown garbage bins mean that batteries must not be thrown into the household waste. The signs below the garbage bins are acronyms for the chemical substances inside the battery:

- Hg: Battery contains Mercury
- Cd: Battery contains Cadmium
- Pb: Battery contains Lead

Technical data*

Sensor	Fuel cell sensor
Accuracy	+/- 0.05‰ (per mill) with reference to 1.00‰ (per mill = default setting in the EU) OR +/- 0.005% BAC with reference to 0.100% BAC (default setting in the USA) OR +/- 0.025mg/L with reference to 0.50mg/L (default setting in UK)
Detection Range	0.00 up to 5.00‰ (per mill) OR 0.000 up to 0.500% BAC (USA) OR 0.00 up to 2.50mg/L (UK)
Warm-up Time	Within 20 seconds at 0.50‰ (EU setting), it may vary depending on the value of per mill which was measured in a test shortly done before
Dimension	120 x 48 x 20 mm (H x W x D)
Weight	93g (batteries included)
Power Supply	2 pieces of 1.5V AAA size Alkaline battery
Working Temperature	5°C - 40°C
Storage Temperature	0°C - 50°C
Calibration	After 500 tests or approx. every 12 months (depending which incident happens first). For a calibration service you can send the device to ACE Instruments (see www.alkomat.net).

* The technical data are subject to change without prior notice for functional improvements.

WARRANTY AND LIABILITY

The ACE Handels- und Entwicklungs GmbH (ACE Instruments) provides a warranty of two (2) years after purchase for material and processing errors for the product (see attached warranty and liability declaration) in all versions and designations; and a warranty of one (1) year for the sensor. This warranty includes the electronics as well as the casing or sensor without sampling components. Components that are inevitably destroyed upon effectiveness (such as fuses) are not part of this warranty. The warranty is valid from the date of delivery or from the date of first activation (always the shorter period), if not stated otherwise in the accompanying ACE Instruments product literature.

MISUSE

OTHER APPLICATIONS OF, OR MODIFICATIONS TO, THIS BREATHALYZER CONSTITUTE A MISUSE AND ENTAIL A RISK OF PERSONAL INJURY AND/OR DAMAGE TO THE DEVICE. ACE INSTRUMENTS SHALL ASSUME NO LIABILITY FOR INJURY/DAMAGE CAUSED BY MISUSE OF THE DEVICE. THE DEVICE IS NOT INTENDED FOR COMMERCIAL USE.

ACE INSTRUMENTS SHALL NOT BE HELD LIABLE FOR ANY DAMAGES THAT OCCUR AS A RESULT OF ANY UNAUTHORIZED MODIFICATIONS, REENGINEERINGS OR AMENDMENTS MADE TO THIS BREATHALYZER OR ITS IMPROPER USE.

LIMITATION OF WARRANTY AND LIABILITY

1. PRODUCT DESCRIPTION

ACE INSTRUMENTS REJECTS ANY FURTHER EXPRESSLY STATED OR IMPLIED WARRANTY, INCLUDING THE MARKETABILITY AND SUITABILITY FOR A DETERMINED PURPOSE. THIS WARRANTY ESPECIALLY EXCLUDES CALIBRATION. AFTER CHANGING THE SENSOR, CALIBRATION IS ABSOLUTELY NECESSARY.

WHEN AN EXCHANGE OF THE SENSOR OR A CALIBRATION IS NECESSARY FOR A PRODUCT, THE USER MAY RECEIVE AN ESTIMATE OF COSTS. SHOULD THE USER REFUSE THE ESTIMATE, ACE INSTRUMENTS WILL SEND THE PRODUCT BACK UNREPAIRED / UNCALIBRATED AND THE PRODUCT WILL BE LABELLED ACCORDINGLY.

SHOULD THE PRODUCT NOT FULFIL THE WARRANTY, THE CLAIM OF THE BUYER AND THE OBLIGATIONS OF ACE INSTRUMENTS ARE LIMITED EXCLUSIVELY AND AT THE DISCRETION OF ACE INSTRUMENTS TO AN EXCHANGE OR REPAIR OF THE AFFECTED GOODS OR TO A REIMBURSEMENT OF THE PURCHASING PRICE FOR THE FAULTY GOODS.

2. APPLICATION RESTRICTION

ACE INSTRUMENTS IS ON NO ACCOUNT LIABLE FOR ADDITIONAL DAMAGES, INCIDENTAL DAMAGES OR CONSEQUENTIAL DAMAGES SUCH AS, AMONG OTHERS BUT NOT LIMITED TO, LOST PROFIT AND DOWNTIME THAT RESULT FROM SELLING, PRODUCING OR USING PRODUCTS, WHICH ARE PART OF THIS APPLICATION RESTRICTION. THIS IS INDEPENDENT OF ANY CLAIMS, SUCH AS FROM CONTRACTS, NEGLIGENCE OR OTHER CULPABLE ACTIONS.

THE PRODUCT, WHICH IS ACCOMPANIED BY THIS WARRANTY AND LIABILITY DECLARATION, NO MATTER IN WHICH VERSION OR DESIGNATION, CAN ESPECIALLY NOT BE USED LEGALLY BINDING TO DETERMINE CAPABILITY TO OPERATE VEHICLES OR TO EXECUTE ACTION THAT COULD LEAD TO DAMAGE TO ONESELF OR OTHERS UNDER THE INFLUENCE OF ALCOHOL.

The thorough check for damage of all products upon receipt, a calibration suitable for the given use by the buyer and the proper use, maintenance and service according to the instructions of the accompanying product literature by ACE Instruments are mandatory prerequisites for this warranty. Maintenance and service work by non-qualified persons, as well as using non-approved consumables or spare parts are a misuse, and all warranty right shall seize immediately. As for all technically challenging products, it is a fundamental prerequisite and mandatory for the warranty of ACE Instruments that all users of this product are familiar with the properties and restrictions of the product as stated in the respective literature.

DISCLAIMER

TO THE EXTENT PERMITTED BY LAW, THIS WARRANTY AND THE REMEDIES SET FORTH ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, REMEDIES AND CONDITIONS, WHETHER ORAL OR WRITTEN, STATUTORY, EXPRESS OR IMPLIED. AS PERMITTED BY APPLICABLE LAW, ACE INSTRUMENTS SPECIFICALLY DISCLAIMS ANY AND ALL STATUTORY OR IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE AND WARRANTIES AGAINST HIDDEN OR LATENT DEFECTS. IF ACE INSTRUMENTS CANNOT LAWFULLY DISCLAIM STATUTORY OR IMPLIED WARRANTIES THEN TO THE EXTENT PERMITTED BY LAW, ALL SUCH WARRANTIES BE LIMITED IN DURATION TO THE DURATION OF THE EXPRESS WARRANTY AND TO THE REPAIR OR REPLACEMENT SERVICE AS STATED ABOVE. EXCEPT AS PROVIDED IN THIS WARRANTY AND TO THE MAXIMUM EXTENT PERMITTED BY LAW, ACE INSTRUMENTS IS NOT RESPONSIBLE FOR DIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ANY LEGAL THEORY.

Some countries, states and provinces do not allow the exclusion or limitation of incidental or consequential damages or allow limitations on how long an implied warranty or condition may last, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary by country, state or province. This Limited Warranty is governed by and construed under the laws of the country in which the product purchase took place.

3. ADVICE RESTRICTION

The buyer confirms that he has checked the intended use and suitability of purchased goods for any purpose in sole responsibility. Between the two parties it is expressly agreed upon that all technical or other advice ACE Instruments regarding use of the goods or services is free of charge and was made with to the exclusive risk of the buyer. Therefor no obligation or liability can arise to ACE Instruments from any recommendation or their result.

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ACE Handels- & Entwicklungs GmbH