ACE AL5500 plus

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Operating Manual

Author's Liability:

All information herein was created conscientiously and under the assumption of correctness. Hence, the author cannot be held liable for the use of any kind whatsoever and/or any ensuing consequences. Furthermore, the author is not liable for any costs or expenses for damage that result from any information in this document. Any information is subject to change without prior notice.

Introduction:

Thank you for purchasing our ACE AL5500 plus Breathalyzer. It measures alcohol concentration in human breath using a semiconductor sensor and calculates the blood alcohol level. The device must only be used pursuant to the instructions mentioned in this manual. In case of any questions you may contact the distributor (ACE) by email to info@ace-instruments.net, www.ace-instruments.net

Note to country-specific measurement units:

Customers in the EU get the breathalyzer with ‰ (per mill) as default setting, in the UK the default setting is mg/L and in the USA it is % BAC. In the chapter "Technical data" you will find information about the accuracy concerning the different measurement units. Please check e.g. the website of the department of transportation, the website of the Ministry of Transport or of automotive clubs to get informed about the current alcohol limits in your country.

Short description:

The ACE AL5500 plus Breathalyzer can be operated by only one button. The test can be done in three different ways: By means of using a mouth piece (= best accuracy), by means of using a straw (= good accuracy) or by just blowing your breath directly upon the sensor opening (= sufficient accuracy). The test results have a tolerance of +/- 0.10‰ (per mill = default setting) with reference to 0.50‰. For devices with UK measurement unit the tolerance is +/- 0.05mg/L with reference to 0.25mg/L, for devices with US measurement unit the tolerance is +/- 0.01% BAC with reference to 0.05% BAC. The test result

is shown on a display with 3 digits (2 decimals) and it can be interpreted 1:1 - so if the shown test result is, e.g., 0.36 then this means 0.36‰ = per mill (for UK: if the shown test result is, e.g., 0.18, then this means 0.18mg/L, for USA: if the shown test result is, e.g., 0.036, then this means 0.036% BAC).


The semiconductor sensor inside the ACE AL5500 plus Breathalyzer can be calibrated. To maintain the measurement accuracy of this breathalyzer it should be calibrated regularly. For a calibration service you can send the device to ACE Instruments (see www.ace-instruments.net).

For each breath test a new mouthpiece is recommended (for hygienic reasons and also to avoid an influence on the test results - a used mouthpiece may contain rest alcohol).

Scope of delivery:

- 1 ACE AL5500 plus Breathalyzer with semiconductor sensor
- 2 AA batteries (1.5 V, Alkaline), already inside the device
- 5 Mouthpieces
- 1 Quick Start Guide ■ 1 Operating Manual (multilingual)

Label of individual components:

This way you receive results that are as accurate as possible:

- A)** Wait at least 15 minutes after your last sip of alcohol. Do not eat, drink or smoke during that time.
- B)** Do not open the casing.
- C)** Do not keep the breathalyzer in cold, hot, moist or badly ventilated places.
- D)** Keep the breathalyzer out of the reach of children.
- E)** Conduct three inlet tests after long storage before conducting the actual alcohol test.
- F)** Do not drop the breathalyzer.
- G)** Change the battery if it is depleted. Diminished battery capacity leads to incorrect test results.
- H)** Do not conduct tests in loud (sound waves disrupt the sensor), drafty or poorly ventilated places.
- I)** After a test result of 1.000‰ (per mill = default setting) or higher the sensor needs rest time. In this case wait at least three minutes before the next test.
- J)** Do not conduct tests in very cold, very hot or an unventilated, smoky ambience.
- K)** Store the breathalyzer in a way to avoid penetration by dust or liquids.
- L)** Blowing too hard may cause damage to the sensor as saliva flows into it.

ACE INSTRUMENTS RECOMMENDS:

- ALWAYS USE A NEW MOUTHPIECE FOR EACH PERSON AND TEST THAT IS CONDUCTED. ACE INSTRUMENTS ASSUMES NO LIABILITY FOR ANY DAMAGES WHICH ARISE WHEN SEVERAL PEOPLE USE THE SAME MOUTHPIECE FOR MULTIPLE TESTS.
- ONLY USE THE PROVIDED SINGLE PACKED AND HYGIENIC MOUTHPIECES - MORE ARE AVAILABLE AT WWW.ACE-INSTRUMENTS.NET.
- DURING EACH TEST THE BREATHALYZER'S INTERNAL AND EXTERNAL PARTS COME INTO CONTACT WITH HUMAN SECRETION (SALIVA, BLOOD, SWEAT, ETC.)! USE DISPOSABLE GLOVES WHEN PULLING OUT A MOUTHPIECE USED BY ANOTHER PERSON! YOU CAN BUY APPROPRIATE GLOVES E.G. AT ACE TECHNIK (WWW.ACE-TECHNIK.DE).

Conducting the measurement:

1. If you use the breathalyzer for the first time, then push the battery compartment cover on the back downwards. Remove the small paper inside the battery compartment (this paper prevents the 2 inserted batteries from discharging). Close the

battery compartment cover again.

2. Insert a mouthpiece (best accuracy) or a straw (good accuracy) into the sensor / mouth opening on the front of the breathalyzer. Make sure that the mouthpiece / straw is firmly seated. Use a new mouthpiece / new straw for each test person.
3. Push the power button*. The device warms up and therefore a countdown starts. Wait until you can see "0" on the display. As soon as you can see "0", start blowing into the mouthpiece.
4. The device will start beeping as soon as you start blowing, keep blowing steadily until you hear a quick double beep tone. Tip: Blow like you would whistle (without a break). If the breath sample was too weak then you can see „Flo“ on the display. Please put off the device and start it again for a new test.
5. After a successful breath sample you will see the digit „0“ moving left to right in the display as the device analyzes the breath sample. The test result is displayed after a few seconds (blinking 10 times). It shows a 3 digit readout (with 2 decimals) of your blood alcohol level in per mill. You can interpret the test result 1:1. The test result has got a tolerance of +/- 0,10 per mill with reference to 0,50 per mill.
6. You can put the device off by pushing the power button until you see "Out" on the display and additionally you will hear a beep tone, but the device also puts itself off automatically approx. 15 seconds after showing the blinking test result.

*If you constantly keep pushing the power button, you can see the number of all tests conducted with the device (starting with zero after a calibration) and also the last result of measurement where alcohol was detected.

IMPORTANT NOTICE:

WHEN YOU HAVE EATEN, DRUNKEN ALCOHOL OR SMOKED A CIGARETTE, WAIT AT LEAST 15 MINUTES BEFORE STARTING THE MEASUREMENT. ANY TEST CONDUCTED PRIOR TO THE EXPIRATION OF THE RECOMMENDED WAITING PERIOD WILL RESULT IN INACCURATE RESULTS.

IF THE TEST RESULT IS BELOW 0.10‰ (PER MILL) OR 0.050 mg/L (UK setting) OR 0.010% BAC (US setting) THEN THE DEVICE SHOWS 0.00 (‰ OR mg/L OR % BAC – depending on the setting) ON THE DISPLAY. THIS HAPPENS BECAUSE AT A LOW RANGE OF ALCOHOL CONCENTRATION IT IS POSSIBLE TO MEASURE OTHER MATERIALS IN YOUR MOUTH.

This breathalyzer is a high quality and sensitive measurement device. To maintain the accuracy of the semiconductor sensor, it should regularly be used. If the device is not used for an extended period of time, the sensor can dry out. ACE Instruments recommends to use the breathalyzer at least 2 times per month (even if you have not been drinking alcohol).

Due to changes in the product, the procedures stated can change. This can take place without prior notice or adaptation of the operating manual and is no cause for complaint, return or price reduction. The user is responsible for recognizing discrepancies between the operating manual and the actual handling of the breathalyzer and must react accordingly.

Only the usage of original ACE mouthpieces

- Guarantees test result accuracy
- Guarantees 100% hygienity
- Guarantees full warranty of this device

ACE INSTRUMENTS ADVISES:

EVEN IF YOU GET A RESULT MEANING „SAFE RANGE“, YOU SHOULD NEVER DRIVE AFTER CONSUMING ALCOHOL OR DO ANY OTHER ACTIVITIES THAT ARE SEVERELY IMPAIRED BY ALCOHOL CONSUMPTION AND MIGHT ENDANGER YOUR HEALTH AND SAFETY OR THE HEALTH AND SAFETY OF OTHERS.

Display indications:

A dark grey rounded rectangle containing the word "Flo" in a white, stylized font.

The breath sample was not sufficient, please repeat the testing process ensuring a steady breath sample (it takes five seconds, just blow like you would whistle without interruption).

A dark grey rounded rectangle containing the word "Cal" in a white, stylized font.

„CAL“ means that the device should be calibrated to maintain accurate measurements.

For a calibration service you can send the device to ACE Instruments (see www.ace-instruments.net).

A dark grey rounded rectangle containing the word "bat" in a white, stylized font.

The battery power is low, please replace the batteries.

A dark grey rounded rectangle containing the word "Out" in a white, stylized font.

The device turns off (either by pressing the power button constantly for 4 seconds or if you don't use the device for about 15 seconds).

Calibration:

To maintain the accuracy of the breathalyzer, ACE recommends a calibration approx. every 6 months. For a calibration service you can send the device to ACE Instruments (see www.ace-instruments.net). The calibration is not mandatory, but the accuracy of the breathalyzer will decrease steadily without it. **We recommend sending it as a parcel, so that in cases of loss an inquiry or proof of delivery is possible.**

Battery notices:

The device associated with this operating manual contains two AA Alkaline batteries (1,5 V) in delivery condition. These batteries are not made for recharging. The use of rechargeable batteries is not possible with this device and can damage it. If the included batteries are depleted, they must not be thrown into the household waste. Used batteries may contain harmful substances that can cause environmental or health damage. Please return all depleted batteries to electrical appliance stores, who collect them for free.

Battery exchange:

To replace the batteries, push the battery cover on the back of the unit downwards and replace the discharged batteries with two new ones. Be aware of the polarity when inserting the batteries (+/- pole). Put the battery compartment cover on the rear panel and press it gently on the device until it clicks. In the EU the return of used batteries is mandatory by law. All depleted batteries will be recycled. This way valuable resources like iron, zinc and nickel can be retrieved and are reusable then. Battery recycling is the easiest lesson in environmental protection - please support it!

The shown garbage bins mean that batteries must not be thrown into the household waste.

The signs below the garbage bins are acronyms for the chemical substances inside the battery:

Hg

Cd

Pb

- Hg: Battery contains Mercury
- Cd: Battery contains Cadmium
- Pb: Battery contains Lead

Technical data*:

Sensor	Semiconductor sensor
Accuracy	+/- 0.10‰ (per mill) with reference to 0.50‰ (per mill = default setting in the EU) OR +/- 10% with reference to 0.05% BAC (default setting in the USA) OR +/- 0.05mg/L with reference to 0.25mg/L (default setting in UK)
Detection Range	0.000 - 4.000‰ (per mill) OR 0.000 - 0.400% BAC (USA) OR 0.000 - 2.000mg/L (UK)
Warm-up Time	Within 30 seconds at 0.5‰ (EU setting), it may vary depending on the value of per mill which was measured in a test shortly done before
Dimension	110 x 46 x 20 mm (H x W x D)
Weight	89g (batteries included)
Power Supply	2 pieces of 1.5V AA size Alkaline battery
Working Temperature	0°C - 40°C
Storage Temperature	0°C - 50°C
Calibration	After approx. every 6 months (Recommendation) For a calibration service you can send the device to ACE Instruments (see www.ace-instruments.net).

***The technical data are subject to change without prior notice for functional improvements.**

WARRANTY AND LIABILITY

WARRANTY

ACE Instruments (ACE) provides a warranty of two (2) years after purchase for material and processing errors for the product (see attached warranty and liability declaration) in all versions and designations; and a warranty of one (1) year for the sensor. This warranty includes the electronics as well as the casing or sensor without sampling components. Components that are inevitably destroyed upon effectiveness (such as fuses) are not part of this warranty. The warranty is valid from the date of delivery or from the date of first activation (always the shorter period), if not stated otherwise in the accompanying ACE product literature.

MISUSE:

OTHER APPLICATIONS OF, OR MODIFICATIONS TO, THIS BREATHALYZER CONSTITUTE A MISUSE AND ENTAIL A RISK OF PERSONAL INJURY AND/OR DAMAGE TO THE DEVICE. ACE SHALL ASSUME NO LIABILITY FOR INJURY/DAMAGE CAUSED BY MISUSE OF THE DEVICE. THE DEVICE IS NOT INTENDED FOR COMMERCIAL USE.

ACE SHALL NOT BE HELD LIABLE FOR ANY DAMAGES THAT OCCUR AS A RESULT OF ANY UNAUTHORIZED MODIFICATIONS, REENGINEERINGS OR AMENDMENTS MADE TO THIS BREATHALYZER OR ITS IMPROPER USE.

LIMITATION OF WARRANTY AND LIABILITY

1. Product description:

ACE REJECTS ANY FURTHER EXPRESSLY STATED OR IMPLIED WARRANTY, INCLUDING THE MARKETABILITY AND SUITABILITY FOR A DETERMINED PURPOSE. THIS WARRANTY ESPECIALLY EXCLUDES CALIBRATION. AFTER CHANGING THE SENSOR, CALIBRATION IS ABSOLUTELY NECESSARY. IF AN EXCHANGE OF THE SENSOR OR A CALIBRATION IS REQUIRED FOR A PRODUCT, THE USER MAY RECEIVE A COST ESTIMATE. SHOULD THE USER REJECT THE ESTIMATE, ACE WILL SEND THE PRODUCT BACK UNREPAIRED/ UNCALIBRATED AND THE PRODUCT WILL BE LABELLED ACCORDINGLY.

SHOULD THE PRODUCT NOT FULFIL THE WARRANTY, THE CLAIM OF THE BUYER AND THE OBLIGATIONS OF ACE ARE LIMITED EXCLUSIVELY AND AT THE DISCRETION OF ACE TO AN EXCHANGE OR REPAIR OF THE AFFECTED GOODS OR TO A REIMBURSEMENT OF THE PURCHASING PRICE FOR THE FAULTY GOODS.

2. Application restriction:

ACE IS ON NO ACCOUNT LIABLE FOR ADDITIONAL DAMAGES, INCIDENTAL DAMAGES OR CONSEQUENTIAL DAMAGES SUCH AS, AMONG OTHERS BUT NOT LIMITED TO, LOST PROFIT AND DOWNTIME THAT RESULT FROM SELLING, PRODUCING OR USING PRODUCTS, WHICH ARE PART OF THIS APPLICATION RESTRICTION. THIS IS INDEPENDENT OF ANY CLAIMS, SUCH AS FROM CONTRACTS, NEGLIGENCE OR OTHER CULPABLE ACTIONS.

THE PRODUCT, WHICH IS ACCOMPANIED BY THIS WARRANTY AND LIABILITY DECLARATION, NO MATTER IN WHICH VERSION OR DESIGNATION, CAN ESPECIALLY NOT BE USED LEGALLY BINDING TO DETERMINE CAPABILITY TO OPERATE VEHICLES OR TO EXECUTE ACTION THAT COULD LEAD TO DAMAGE TO ONESSELF OR OTHERS UNDER THE INFLUENCE OF ALCOHOL.

THE THOROUGH CHECK FOR DAMAGE OF ALL PRODUCTS UPON RECEIPT, A CALIBRATION SUITABLE FOR THE GIVEN USE BY THE BUYER AND THE PROPER USE, MAINTENANCE AND SERVICE ACCORDING TO THE INSTRUCTIONS OF THE ACCOMPANYING PRODUCT LITERATURE BY ACE ARE MANDATORY PREREQUISITES FOR THIS WARRANTY. MAINTENANCE AND SERVICE WORK BY NON-QUALIFIED PERSONS, AS WELL AS USING NON-APPROVED CONSUMABLES OR SPARE PARTS ARE A MISUSE, AND ALL WARRANTY RIGHT SHALL SEIZE IMMEDIATELY. AS FOR ALL TECHNICALLY CHALLENGING PRODUCTS, IT IS A FUNDAMENTAL PREREQUISITE AND MANDATORY FOR THE WARRANTY BY ACE THAT ALL USERS OF THIS PRODUCT ARE FAMILIAR WITH THE PROPERTIES AND RESTRICTIONS OF THE PRODUCT AS STATED IN THE RESPECTIVE LITERATURE.

DISCLAIMER

TO THE EXTENT PERMITTED BY LAW, THIS WARRANTY AND THE REMEDIES SET FORTH ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, REMEDIES AND CONDITIONS, WHETHER ORAL OR WRITTEN, STATUTORY, EXPRESS OR IMPLIED. AS PERMITTED BY APPLICABLE LAW, ACE SPECIFICALLY DISCLAIMS ANY AND ALL STATUTORY OR IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE AND WARRANTIES AGAINST HIDDEN OR LATENT DEFECTS. IF ACE CANNOT LAWFULLY DISCLAIM STATUTORY OR IMPLIED WARRANTIES THEN TO THE EXTENT PERMITTED BY LAW, ALL SUCH WARRANTIES BE LIMITED IN DURATION TO THE DURATION OF THE EXPRESS WARRANTY AND TO THE REPAIR OR REPLACEMENT SERVICE AS STATED ABOVE. EXCEPT AS PROVIDED IN THIS WARRANTY AND TO THE MAXIMUM EXTENT PERMITTED BY LAW, ACE IS NOT RESPONSIBLE FOR DIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ANY LEGAL THEORY.

Some countries, states and provinces do not allow the exclusion or limitation of incidental or consequential damages or allow limitations on how long an implied warranty or condition may last, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary by country, state or province. This Limited Warranty is governed by and construed under the laws of the country in which the product purchase took place.

3. Advice restriction:

The buyer confirms that he has checked the intended use and suitability of purchased goods for any purpose in sole responsibility. Between the two parties it is expressly agreed upon that all technical or other advice by ACE regarding use of the goods or services is free of charge and was made with to the exclusive risk of the buyer. Therefore no obligation or liability can arise to ACE from any recommendation or their result.

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La empresa ACE Handels- und Entwicklungs GmbH en ningún
ACE Handels und Entwicklungs GmbH